

## 5'-Nucleotidase Assay Kit (5'-NT)

**Method:** Colorimetric (Kinetic)

Cat .No.	Size	Instrument
GB8010G	R1: 2 × 20 ml R2: 1 × 20 ml	For Hitachi 717 & ShimadzuCL7200/8000
GS8011G	R1: 2 × 20 ml R2: 1 × 20 ml	For Hitachi917 & OlympusAU640/400/600
GX8011G	R1: 2 × 20 ml R2: 1 × 20 ml	For SYNCHRON CX4-5-7-9/LX20

### INTENDED USE



For the *in vitro* quantitative determination of 5'-nucleotidase activity in human serum or plasma.

### CLINICAL SIGNIFICANCE

5'-NT is an enzyme catalyzing the hydrolysis of nucleoside-5'-monophosphates to nucleosides and inorganic phosphate. The enzyme is widely distributed in human and animal tissues. The activity present in sera is released from the membrane of liver cells by bile salts and has been used as a marker for liver disease<sup>[1]</sup>. Increased enzyme levels in sera are associated with certain forms of liver disease, such as intra or extra-hepatic obstruction and particularly in cases of hepatic carcinoma as well as in mastectomy patients with recurrent metastases. The diagnostic value of 5'-NT has been shown to be superior to other liver enzymes, especially in cases of liver metastasis.

### ASSAY PRINCIPLE

The 5'-NT assay is based on the enzymatic hydrolysis of 5'-monophosphate (5'-IMP) to form inosine which is converted to hypoxanthine by purine nucleoside phosphorylase (PNP). Hypoxanthine is then converted to uric acid and hydrogen peroxide (H<sub>2</sub>O<sub>2</sub>) by xanthine oxidase (XOD). H<sub>2</sub>O<sub>2</sub> is further reacted with N-Ethyl-N-(2-hydroxy-3-sulfopropyl)-3-methylaniline (EHSPT) and 4-aminoantipyrine (4-AA) in the presence of peroxidase (POD) to generate quinone dye which is monitored in a kinetic manner. The entire enzymatic reaction scheme is shown below.



### SAMPLE COLLECTION AND PREPARATION

Serum or plasma samples. Use fresh and non-hemolyzed serum or EDTA treated plasma samples. Serum or plasma samples are stable for a week at 4°C.

### REAGENT COMPOSITION

Contents	Concentration of Solutions
<b>Reagent 1 (R1)</b>	
Goods buffer (PH7.6)	100 mmol/L
4-AA	2 mmol/L
PNP	0.1 U/L
XOD	0.2 U/L
POD	0.6 U/L
Stabilizer	
<b>Reagent 2 (R2)</b>	
Goods buffer (PH7.6)	100 mmol/L
IMP	10 mmol/L
EHSPT	2 mmol/L

### STABILITY AND PREPARATION OF REAGENTS

All reagents are ready to use. Stable up to the expiry date when stored at 2-8°C. The assay kit reagents are stable for 1 month after opening and kept at 2-8°C.


### ASSAY PROCEDURE

**Test Procedure for Analyzers** (HITACHI 7170/917)

Assay Mode: 2 Point Rate 28-34

Wave Length (main/sub): 546 nm/800 nm

Sample: 5 μl  
R1: 180 μl R2: 90 μl



- Mix 5 μl sample with 180 μl R1 and incubate at 37°C for 5 minutes.
- Add 90 μl R2 into cuvette, mix and incubate for 3 minutes at 37°C.
- Read initial absorbance and start timer simultaneously, read again after 1 and 2 minutes.
- Calculate absorbance change per minute (ΔA/min)

## CALCULATION ( $\varepsilon=16.18$ )

$$5'-NT \text{ (U/L)} = \frac{\Delta A/\text{min} \times V_t}{\varepsilon \times V_s \times L} = \Delta A/\text{min} \times 3400$$

## QUALITY CONTROL

For quality control, use 5'-NT control as daily quality control sera and can be purchased separately. Values should fall within a specific range. If these values fall outside the range and repetition excludes error, the following steps should be taken:

1. Check instrument settings and light source.
2. Check reaction temperature.
3. Check expiration date of kit and contents.

## REFERENCE VALUE

Serum or plasma: 0-10 U/L

It is recommended that each laboratory establish its own reference range to reflect the age, sex, diet and geographical location of the population.

## LINEARITY

The method is linear up to 300 U/L. Samples above this concentration should be diluted 1+1 with 0.9% NaCl and reassayed. Multiply the result by 2.

## SPECIFIC PERFORMANCE CHARACTERISTICS INTERFERENCE

Hemoglobins  $\leq 500$  mg/dl, Intralipid  $\leq 1000$  mg/dl  
ALP  $\leq 1250$  U/l, without interference.  
Bilirubin  $\leq 40$  mg/dl, Ascorbic Acid  $\leq 50$  mg/dl, interference.

## PRECISION

The CV of the test should be  $\leq 10\%$ .

Intra assay precision		
N=20	Level1	Level 2
Mean	7.5	78.8
SD	0.11	0.95
CV	1.5%	1.2%
Inter assay precision		
N=5	Level1	Level 2
Mean	7.83	79.3
SD	0.16	1.35
CV	2.0%	1.7%

## SENSITIVITY

The minimum detectable level that can be distinguished from zero has been determined as 1.81 U/L.

## CORRELATION

This method (Y) was compared with another commercially available method (X) and the following linear regression equation obtained:  
 $Y=1.0426X+1.3074$ ,  $R^2=0.9912$ ; 64 patient samples were analyzed.

## SAFETY PRECAUTIONS AND WARNINGS

1. For *in vitro* diagnostic use only. Do not pipette by mouth. Exercise the normal precautions required for handling laboratory reagents.
2. Reagents contain Sodium Azide. Avoid ingestion or contact with skin or mucous membranes. In case of skin contact, flush affected area with copious amounts of water. In case of contact with eyes or if ingested, seek immediate medical attention.
3. Sodium Azide reacts with lead and copper plumbing, to form potentially explosive azides. When disposing of such reagents flush with large volumes of water to prevent azide build up. Exposed metal surfaces should be cleaned with 10% sodium hydroxide.
4. All specimens used in this test should be considered potentially infectious. Universal Precautions, as they apply at your facility, should be used for handling and disposing of materials during and after testing.

## REFERENCES

1. Eroglu A. Activities of adenosine deaminase and 5' - nucleotidase in cancerous and noncancerous human colorectal tissues[J]. Med Oncol 2000; 17(4):319 - 24.1
2. Alain Bertrand and Jean Buret, A one-step determination of 5'-nucleotidase using a centrifugal analyzer. Clinica Chimica Acta, 119(1982)275-284.
3. By Z. Ahmed and J.L.Reis, The Activation and Inhibition of 5-Nucleotidase. Clin Chem. 1998, 69 (11), 102-106.

## INDEX OF SYMBOLS



Manufacture



Catalogue Number

Lot number



Date of manufacture



Use by (Expiration date)



For In-Vitro Diagnostic use only



Stored at 2-8°C



Attention: See instruction for use



Authorized Representative in the

European Company

## INSTRUMENT SETTINGS FOR HITACHI 917

Hitachi 7170 Parameter Application

Gcell

5'-NT  
Cat. No: GB8010G/GS8011G/GX8011G

<b>Analysis</b>		<b>Range</b>	
Test / Type	SNT	Application Code	SNT
Assay/Time/Point	Rate A	Report Name	SNT
Wave (Sub/Main)	800 A 546 A	Data Mode	On Board A
S.Vol (Normal)	5 0.0 0	Control Interval	0
S.Vol (Decrease)	2 0.0 0	Instrument Factor (Y=aX+b)	a= 1.0 b= 0
S.Vol (Increase)	10 0.0 0	Technical Limit	0 300
Diluent	Water 0	Expected Value	
Reagent (R1) T1	180 0 *	Qualitative	
Reagent (R2) T2	0 0 00000 0	(Male)	0 Y A
Reagent (R3) T3	90 0 *	(Female)	0 Y A
Reagent (R4) T4	0 0 00000 0		
Abs. Limit	32000 Increase A		
Prozone Limit	0 34 Lower A		
Cell Detergent	Detergent 1 A		

<b>Calibration</b>		<b>STD Conc</b>	
Calibration type	Linear A A	<Standard>	(1) (2) (3) (4) (5) (6)
Point	1 Span Point 0	Concentration	0.0 0 0 0 0 0
Weight	0	Position	Water 0 0 0 0 0 0
Auto calibration		Volume	5 0 0 0 0 0
Time Out		<Pre-Diluent>	
Blank	0	Volume	0 0 0 0 0 0
Span	0	Diluent	0 0 0 0 0 0
2Point	0	Cal. Code	0 0 0 0 0 0
Full	0		
SD Limit	999		
Duplicate limit	1900		
Sensitivity limit	0		
S1 Abs limit	-32000 32000		

Attention: \* entered by operator  
K-factor = 34000

## INSTRUMENT SETTINGS FOR HITACHI 902

Hitachi 7020 Instrument Settings

Gcell

5'-NT  
Cat. No: GB8010G/GS8011G/GX8011G

No.	<Chemistry>	46	Sens. Limit	0
1	Test Name	47	S 1 ABS Limit (L)	-32000
2	Assay Code (Mthd)	48	S 1 ABS Limit (H)	32000
3	Assay Code (2. Test)	49	ABS Limit	32000
4	Reaction Time	50	ABS Limit (D/I)	Increase
5	Assay Point 1	51	Prz. Limit	0
6	Assay Point 2	52	Prz. Limit (U/D)	Lower
7	Assay Point 3	53	Prz. (End Point)	35
8	Assay Point 4	54	Expect. Value (L)	0
9	Wave Leng. (SUB)	55	Expect. Value (H)	10
10	Wave Leng. (MAIN)	56	Instr. Fact. (a)	1
11	Sample Volume	57	Instr. Fact. (b)	0
12	R1 VOLUME	58	Key Setting	*
13	R1 Pos.			
14	R1 Bottle Size			
15	R2 VOLUME			
16	R2 Pos.			
17	R2 Bottle Size			
18	R3 VOLUME			
19	R3 Pos.			
20	R3 Bottle Size			
21	Calib. Type (Type)			
22	Calib. Type (Wght)			
23	Calib. Conc. 1			
24	Calib. Pos. 1			
25	Calib. Conc. 2			
26	Calib. Pos. 2			
27	Calib. Conc. 3			
28	Calib. Pos. 3			
29	Calib. Conc. 4			
30	Calib. Pos. 4			
31	Calib. Conc. 5			
32	Calib. Pos. 5			
33	Calib. Conc. 6			
34	Calib. Pos. 6			
35	S 1 ABS.			
36	K Factor			
37	K 2 Factor			
38	K 3 Factor			
39	K 4 Factor			
40	K 5 Factor			
41	A Factor			
42	B Factor			
43	C Factor			
44	SD Limit			
45	Duplicate Limit			

\* Data entry by the user

## INSTRUMENT SETTINGS FOR Olympus400/640/2700

Olympus AU640/400/2700 Instrument Settings

Gcell

5'-NT  
Cat. No:GB8010G/GS8011G/GX8011G

<b>Specific Test Parameters</b>		<b>Calibration Specific</b>	
Test Name:	SNT	Test No.:	MB
Sample: Volume	5	Cal. Type:	MB
Reagents: R1 Volume	180	Formula:	Y = AX + B
R2 Volume	90	Calibration Selection:	
Dilution	0	Point 1	Cal. No.
Dilution	0	Point 2	OD
Wavelength: Pri.	840	Point 3	Conc.
Method:	RATE	Point 4	Factor/OD-L
Reaction Slope:	-	Point 5	Factor/OD-H
Measuring Point 1: First	27	Point 6	
Measuring Point 2: First	27	Point 7	
Linearity:	30%	1-Point Cal. Point:	
No-Lag-Time:	YES	MB Type Factor:	3400
		Cal. Stability Period:	

Attention: \* Entered By Operator

## INSTRUMENT SETTINGS FOR CX4/5/7/9

Synchron CX-4/5/7/9 User-defined Chemistries

Gcell

5'-NT  
Cat. No: GS8011G/GB8010G/GX8010G

<b>USER ID:</b>		<b>Calculate Factor:</b> 6800	
<b>Chemistry Name:</b> 5'-NT		<b>Math Model:</b> Linear	
<b>Test Name:</b> 5'-NT		<b>Cal Time Limit:</b> 0 Hrs	
<b>Reaction Type:</b> Rate 1		<b>No. Of Calibrators:</b> 0	
<b>Reaction Direction:</b> Positive			
<b>Units:</b> U/L			
<b>Decimal Precision:</b> X.X			
<b>Primary Wavelength:</b> 560 nm		<b>Secondary Wavelength:</b> 700 nm	
<b>Sample Volume:</b> 4 µl		<b>CALIBRATORS</b>	
<b>Primary Inject Rgt:</b>		<b>MULTIPOINT SPAN</b>	
A: 144 µl		#1: 1 - 2 0.000	
B: 72 µl		#2: 2 - 3 0.000	
<b>Secondary Inject Rgt:</b>		#3: 3 - 4 0.000	
None: 0 µl		#4: 4 - 5 0.000	
<b>Add Time:</b> sec		#5: 5 - 1 0.000	
<b>RAGENT BLANK</b>		<b>REACTION</b>	
Start Read: 260 sec		Start Read: 180 sec	
End Read: 300 sec		End Read: 300 sec	
Low ABS Limit: -1.500		Low ABS Limit: -1.500	
High ABS Limit: 1.500		High ABS Limit: 1.500	
<b>USABLE RANGE</b>		<b>SUBSTRATE DEPLETION</b>	
Lower Limit: 0		Initial Rate: 99.999	
Upper Limit: 300		Delta ABS: 1.5	
<b>RECOVERY/SENSITIVITY</b>			
Std Dev (conc): *			
CV (%): *			
Std Dev (mA): *			
Threshold: *			

Attention: \* Entered By Operator