

Lactate Dehydrogenase –L Assay Kit (LDH)

Method: DGKC

Cat . No	Size	Instrument
GB230X	R1: 2×80 ml R2: 1×40 ml	For Hitachi 717 &ShimadzuCL7200/8000
GS231X	R1: 2×70 ml R2: 1×35 ml	For Hitachi 917 &OlympusAU640/400/600
GH231X	R1: 2×50 ml R2: 1×25 ml	For Hitachi 902
GX231X	R1: 2×80 ml R2: 2×20 ml	For SYNCHRON CX4-5-7-9 /LX20/DXC600-800
GT231X	R1: 5×48 ml R2: 2×30 ml	For TOSHIBA


INTENDED USE

For the *in vitro* quantitative determination of lactate dehydrogenase activity in human serum.

CLINICAL SIGNIFICANCE^[1,2]

The enzyme lactate dehydrogenase (LDH-L) is distributed in tissues particularly heart, liver, muscle, and kidney. The enzyme found in circulation is a mixture of five isoenzymes based on their mobility. Elevated serum levels of LDH-L are found in serum in myocardial infarction, liver disease, renal disease, certain forms of anemia, malignant diseases and progressive muscle dystrophy.

ASSAY PRINCIPLE^[3]



LDH catalyzes the oxidation of lactate to pyruvate, and NAD is reduced to NADH, which can be measured at 340 nm.

SAMPLE COLLECTION AND PREPARATION

Serum samples.

REAGENT COMPOSITON

Contents	Concentration of Solutions
Reagent 1 (R1)	
N-methyl-D-glucamine	325 mmol/l
Lactate	50 mmol/l
Reagent 2 (R2)	
NAD ⁺	10 mmol/l

STABILITY AND PREPARATION OF REAGENTS

All reagents are ready to use.


Stable up to the expiry date when stored at 2-8°C. The LDH assay kit reagents are stable for 1 month on-board the analyzer after opening and kept at 2-8°C.

ASSAY PROCEDURE

Test Procedure for Analyzers (HITACHI7170/917)

Assay Mode: 2 Point Rate 20-29

Wave Length (main/sub): 340 nm/405 nm



1. Mix 8 µl sample with 200 µl R1 and incubate at 37°C for 5 minutes.
2. Add 50 µl R2 into cuvette, mix and incubate for 1 min at 37°C.
3. Read initial absorbance and start timer simultaneously, read again after 1 and 2 minutes.
4. Calculate absorbance change per minute ($\Delta A/\text{min}$)

CALCULATION

$$\text{Concentration} = \frac{A_{\text{sample}} / \text{min}}{A_{\text{calibrator}} / \text{min}} \times \text{Calibrator value}$$

$$\text{LDH (U/L)} = \frac{\Delta A / \text{min} \times T_v}{\epsilon \times S_v \times L} = \Delta A / \text{min} \times 10080$$

($\epsilon = 6.22$)

CALIBRATION

Recommend that this assay should be calibrated using Randox Calibration Serum Level 3 or Level 2.

QUALITY CONTROL

Randox Assayed Multisera, Level 2 and Level 3 are recommended for daily quality control. Two levels of controls should be assayed at least once a day. Values obtained should fall within a specified range. If these values fall outside the range and repetition excludes error, the following steps should be taken:

1. Check instrument settings and light source.
2. Check reaction temperature.
3. Check expiration date of kit and contents.

NORMAL VALUE

	Female (U/L)	Male (U/L)
Adult	135 - 215	135 - 225
1 - 3 years old	165 - 395	155 - 345
4 - 6 years old	135 - 345	155 - 345
7 - 9 years old	140 - 280	145 - 300
10 -12 years old	120 - 260	120 - 325
13 -15 years old	100 - 275	120 - 290
16-18 years old	105 - 230	105 - 235

It is recommended that each laboratory establish its own reference range to reflect the age, sex, diet and geographical location of the population.

LINEARITY

The method is linear up to 1000 U/L. If the sample above this concentration should be diluted 1+1 with 0.9% NaCL and repeat assay. Multiply the result by 2.

SPECIFIC PERFORMANCE CHARACTERISTICS INTERFERENCE

The following analytes were tested up to the levels indicated and found not to interfere:

Hemoglobin:	50 mg/dl
Intralipid:	1000 mg/dl
Bilirubin:	80 mg/dl
Ascorbic acid:	50 mg/dl

PRECISION

The CV of the test should be less than 5%

Inter assay precision		
N=5	Level1	Level 2
Mean (U/l)	210.4	337.3
SD	2.6	3.1
CV	1.3%	0.92%
Intra assay precision		
N=20	Level1	Level 2
Mean (U/l)	185.3	318.4
SD	1.1	2.3
CV	0.6%	0.7%

SENSITIVITY

The minimum detectable level that can be distinguished from zero has been determined as 5.02 U/L.

CORRELATION

This method (Y) was compared with another commercially available method (X) and the following linear regression equation obtained:
 $Y=0.940X+3.993$, $R^2=0.997$; 60 patient samples were analyzed .

SAFETY PRECAUTIONS AND WARNINGS

1. For in vitro diagnostic use only. Do not pipette by mouth. Exercise the normal precautions required for handling laboratory reagents.
2. The reagents contains Sodium Azide. Avoid ingestion or contact with skin or mucous membranes. In case of skin contact, flush affected area with copious amounts of water. In case of contact with eyes or if ingested, seek immediate medical attention.
3. Sodium Azide reacts with lead and copper plumbing, to form potentially explosive azides. When disposing of such reagents flush with large volumes of water to prevent azide build up. Exposed metal surfaces should be cleaned with 10% sodium hydroxide.
4. All specimens used in this test should be considered potentially infectious. Universal Precautions, as they apply at your facility, should be used for handling and disposing of materials during and after testing.

REFERENCES

1. Amador E., Dorfman L.E., Wacker W.E., Clin. Chem.,1963; 9: 331.
2. Henry, R.J., Clinical Chemistry, Principles and Techniques, 2nd Edition, Harper and Row, p. 819, 1974.
3. Tietz, N., Fundamentals of Clinical Chemistry, W.B.Saunders Co., Philadelphia, PA, p. 652, 1976.

INDEX OF SYMBOLS



Manufacture



Catalogue Number
Lot number



Date of manufacture



Use by (Expiration date)




For In-Vitro Diagnostic use only



Stored at 2-8°C



Attention: See instruction for use



Authorized Representative in the
European Company

INSTRUMENT SETTINGS FOR HITACHI 917

Hitachi 7170 Parameter Application

Gcell

LDH
Cat. No: GS231X-GB230X

Analysis	
Test / Type	LDH
Assay/Time/Point	Rate A A 10 A 20 29 0 0
Wave (Sub/Main)	405 A 340 A
S.Vol (Normal)	8 0.0 0
S.Vol (Decrease)	2 0.0 0
S.Vol (Increase)	8 0.0 0
Diluent	Water 0
Reagent (R1) T1	200 0 * 0
Reagent (R2) T2	0 0 00000 0
Reagent (R3) T3	50 0 * 0
Reagent (R4) T4	0 0 00000 0
Abs. Limit	32000 Increase A
Prozone Limit	0 34 Lower A
Cell Detergent	Detergent 1 A

Range	
Application Code	* Unit U/L A
Report Name	LDH
Data Mode	On Board A
Control Interval	0
Instrument Factor (Y=aX+b)	a=1.0 b=0
Technical Limit	0 1500
Expected Value	
(Male)	0 Y A (1)0
	100 Y A (2)0
	135 225 (3)0
(Female)	0 Y A (4)0
	100 Y A (5)0
	135 215

Calibration	
Calibration type	Linear A A
Point	1 Span Point 0
Weight	0
Auto calibration	
Time Out	
Blank	0 Change Over Blank A
Span	0 Blank A
2Point	0
Full	0
SD Limit	0.1
Duplicate limit	1000
Sensitivity limit	0
S1 Abs limit	-32000 32000

STD Conc	
<Standard>	(1) (2) (3) (4) (5) (6)
Concentration	0 0 0 0 0 0
Position	Water 0 0 0 0 0 0
Volume	8 0 0 0 0 0
<Pre-Diluent>	
Volume	0 0 0 0 0 0
Diluent	0 0 0 0 0 0
Cal. Code	0 0 0 0 0 0

Attention: * entered by operator
K-factor= 5185

INSTRUMENT SETTINGS FOR Olympus400/640/2700

Olympus AU640/400/2700 Instrument Settings

Gcell

LDH
Cat. No: GS231X-GB230X

Specific Test Parameters	
Test Name:	LDH
Sample Volume:	4
Dilution:	0
Pre-Dilution Rate:	0
Reagents: R1 Volume:	200
Dilution:	0
Min OD:	-2.0000
Max OD:	2.5000
Reagent OD Limit:	L -2.0000 H 2.5000
Wavelength: Pri:	340
Sec:	410
Method:	RATE
Reaction Slope:	+
Measuring Point 1: First	14 Last 19
Measuring Point 2: First	19 Last
Linearity:	30%
No-Lag-Time:	YES
Dynamic Range:	L 0 H 1512
Correlation Factor:	A 1.0 B 0.0
Onboard Stability Period:	

Calibration Specific	
Test No.:	MB
Name:	LDH
Type:	SER
Cal. Type:	2
Formula:	Y = AX + B
Process:	CONC
Calibration Selection:	
Cal. No.	OD
Conc.	Factor/OD-L
Factor/OD-H	
Point 1	
Point 2	
Point 3	
Point 4	
Point 5	
Point 6	
Point 7	
1-Point Cal. Point:	
MB Type Factor:	10080
Cal. Stability Period:	

Attention: * Entered By Operator

INSTRUMENT SETTINGS FOR HITACHI 902

Hitachi 7020 Instrument Settings

Gcell

LDH
Cat. No: GB 230X/GS 231X/GH 231X

No.	<Chemistry>
1	Test Name LDH
2	Assay Code (Mthd) Rate A
3	Assay Code (2. Test) 0
4	Reaction Time 10
5	Assay Point 1 20
6	Assay Point 2 30
7	Assay Point 3 0
8	Assay Point 4 0
9	Wave Leng. (SUB) 405
10	Wave Leng. (MAIN) 340
11	Sample Volume 4
12	R1 VOLUME 200
13	R1 Pos. *
14	R1 Bottle Size Large
15	R2 VOLUME 0
16	R2 Pos. 0
17	R2 Bottle Size Small
18	R3 VOLUME 50
19	R3 Pos. *
20	R3 Bottle Size Small
21	Calib. Type (Type) K Factor
22	Calib. Type (Weight) 0
23	Calib. Conc. 1 0
24	Calib. Pos. 1 99
25	Calib. Conc. 2 0
26	Calib. Pos. 2 0
27	Calib. Conc. 3 0
28	Calib. Pos. 3 0
29	Calib. Conc. 4 0
30	Calib. Pos. 4 0
31	Calib. Conc. 5 0
32	Calib. Pos. 5 0
33	Calib. Conc. 6 0
34	Calib. Pos. 6 0
35	S1 ABS. 0
36	K Factor 10080
37	K 2 Factor 10000
38	K 3 Factor 10000
39	K 4 Factor 10000
40	K 5 Factor 10000
41	A Factor 0
42	B Factor 0
43	C Factor 0
44	SD Limit 999
45	Duplicate Limit 1000
46	Sens. Limit 0
47	S1 ABS Limit (L) -32000
48	S1 ABS Limit (H) 32000
49	ABS Limit 32000
50	ABS Limit (D/I) Increase
51	Prz. Limit 0

52	Prz. Limit (U/D) Lower
53	Prz. (End Point) 35
54	Expect. Value (L) 135
55	Expect. Value (H) 225
56	Instr. Fact. (a) 1
57	Instr. Fact. (b) 0
58	Key Setting *

* Data entry by the user

INSTRUMENT SETTINGS FOR CX4/5/7/9

Synchron CX-4/5/7/9 User-defined Chemistries

Gcell

LDH
Cat. No: GB230X/GS231X/GX231X

USER ID:	
Chemistry Name:	LDH
Test Name:	LDH
Calculate Factor:	20160
Reaction Type:	Rate 1
Reaction Direction:	Positive
Units:	U/L
Math Model:	Linear
Cal Time Limit:	336 Hrs
No. of Calibrators:	0
Decimal Precision:	X
Primary Wavelength:	340 nm
Secondary Wavelength:	410 nm
Sample Volume:	4 µl
Primary Inject Rgt:	A: 200 µl
None:	µl
Secondary Inject Rgt:	B: 50 µl
Add Time:	512 sec
CALIBRATORS	MULTIPOINT SPAN
#1:	1 - 2 0.000
#2:	2 - 3 0.000
#3:	3 - 4 0.000
#4:	4 - 5 0.000
#5:	5 - 1 0.000
RAGENT BLANK	
Start Read:	464 sec
End Read:	480 sec
Low ABS Limit:	-1.500
High ABS Limit:	1.500
USABLE RANGE	
Lower Limit:	0
Upper Limit:	1500
RECOVERY/SENSITIVITY	
Std Dev (conc):	*
CV (%):	*
Std Dev (mM):	*
Threshold:	*
REACTION	
Start Read:	256 sec
End Read:	346 sec
Low ABS Limit:	-1.500
High ABS Limit:	1.500
SUBSTRATE DEPLETION	
Initial Rate:	99.999
Delta ABS:	1.5

Attention: * Entered By Operator